

# FCoV Ag

## General Information

**Purpose :** Detection of Feline Corona virus antigen

**Principle :** Immunochromatographic assay

**Specimen :** Feline feces

**Materials provided :**

- Anigen Rapid FCoV test device
- Assay diluents tube
- Disposable swab
- Disposable dropper

**Special Features :**


- One-step testing procedure : Quick & Accurate
- High sensitivity and specificity
  - Sensitivity : 95.5 % vs. PCR
  - Specificity : 98.8 % vs. PCR




◆ **Feline infectious peritonitis ; FIP**

There are two basic forms of FIP-effusive (wet) and non effusive (dry) that have been characterized. Cats with effusive FIP usually only survive a few days-weeks at best. Cats with noneffusive FIP can survive for weeks or months, although after neurologic signs begin death usually occurs rapidly.

## Test Procedure



### Interpretation



## Ordering Information

Cat. No.	Description	Type	Pack size
RG1205DD	FCoV Ag	Device	10 Tests/Kit