

# FPV Ag

## General Information

**Purpose :** Detection of Feline Parvo virus antigen

**Principle :** Immunochromatographic assay

**Specimen :** Feline feces

**Materials provided :**


- Anigen Rapid FPV test device
- Assay diluents tube
- Disposable swab
- Disposable dropper

**Special Features :**

- One-step testing procedure : Quick & Accurate
- High sensitivity and specificity
  - Sensitivity : 97% vs. RT-PCR
  - Specificity : 98.5% vs. RT-PCR



## Test Procedure



## Ordering Information

Cat. No.	Description	Type	Pack size
RG1203DD	FPV Ag	Device	10 Tests/Kit