

# Toxoplasma Ab

## General Information

**Purpose :** Detection of *Toxoplasma gondii* antibody

**Principle :** Immunochromatographic assay

**Specimen :** Whole blood, plasma or serum

**Materials provided :**


- Anigen Rapid Toxo test device
- Assay diluents dropping bottle
- Anticoagulant tube
- Disposable capillary tube

**Special Features :**

- One-step testing procedure : Quick & Accurate
- High sensitivity and specificity
  - Sensitivity : 99.9% vs. ELISA
  - Specificity : 98.7% vs. ELISA



## Test Procedure



## Ordering Information

Cat. No.	Description	Type	Pack size
RB22-06	Toxoplasma Ab	Device	10 Tests/Kit